nithya

Nithya Collagen & Nithya S-Line

Endless
Beauty

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
nithya

Nithya Collagen & Nithya S-Line

Anti-Age Activity and Tolerance Evaluation of Collagen Micro-Injection Treatment Associated to Topical Application of a Cosmetic Formulation (Investigator-Initiated Multicentre Trial)

Adele Sparavigna¹, Antonello Tateo, Elena Inselvini, Massimiliano Tocchio, Maria Cristina Orlandini and Giorgio Botali

¹Dermlng srl, Clinical Research and Bioengineering Institute, Milan, Italy

*Corresponding author: Adele Sparavigna, Dermlng srl, Clinical Research and Bioengineering Institute, Milan, Italy, Tel: +39 3420399117; E-mail: adele.sparavigna@dermlng.com

Received date: March 20, 2017; Accepted date: April 11, 2017; Published date: April 14, 2017

Copyright: © 2017 Sparavigna A. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Abstract

Objective: A novel equine type I collagen treatment consisting of an intradermal injection associated with the use of a topical mousse was developed with characteristics especially suited for global face rejuvenation. A multicenter, investigator-initiated, open clinical trial, conducted by 6 Italian centers, was carried out to evaluate the anti-age performance of this combined treatment.

Methods: The study was conducted on 72 female healthy subjects, age range 40-65 years, asking for midface volume restoration and presenting cutaneous aging/photoaging signs. The 1st injection treatment with the injectable product was performed immediately after the basal evaluation (T0). Two touch-up treatments were performed after 2 (T1) and 4 weeks (T2). Collagen mousse was to be applied by the volunteers on the face (including the submental area) twice a day, with a mild massage. Subject returned 2 (T3), 3 (T4) and 6 months (T5) after the 1st injection treatment for the follow-up evaluation. The aesthetics results were established through the use of clinical evaluations (Wrinkle Severity Rating Scale (WSRS), Facial Volume Loss Scale (FVLS), and wrinkles grade of the area around the eyes (Glogau's reference photographic scale)) and instrumental assessments (profilometry on skin replicas) supported by photographic documentation and face volume image analysis. Secondary endpoints were tolerance evaluation, performance duration and subjects' efficacy judgement.

Results: The study treatment determined a very significant reduction of all clinical parameters considered (Crows' feet grade, FVLS and WSRS) at every study time. Profilometry on crows' feet skin replicas and face volume image analysis confirmed the clinical evaluation, showing a statistical/clinical significant reduction of all the profilometric parameters (anti-wrinkle efficacy) and a consistent improvement of cheek volume (bio-revolumetric effect).

Conclusion: Obtained results confirm the anti-aging activity of the associated collagen treatment (injectable dermal filler and topical mousse). The aesthetic performance resulted persistent up to the final follow-up visit, sign of a long-lasting stimulating activity on cellular functionality of the associated treatment. The majority of volunteers underlined a very marked reshaping of face contour as well as an important improvement of skin smoothness, brightness and hydration. The study treatment was well tolerated, no unexpected adverse reaction related to the tested products/injection procedure occurred during the trial.

Keywords: Equine collagen; Dermal filler; Topical treatment; Bio restructuring; Facial rejuvenation; Volume loss; Wrinkles

Introduction

Background

The physiological aging process, adversely affected by ultraviolet radiation, environmental factors, unbalanced diet, pollution and hormonal factors, causes major changes to the skin [1,2]. The loss of the normal level of epidermis hydration is followed by major structural alterations of the extracellular matrix (ECM), such as degradation of collagen [1-3]. The result is severe structural and functional impairment of the dermis, with inevitable aesthetic consequences: sagging skin, loss of volume and deep wrinkles [4]. Collagen, a fibrous protein of the ECM, is the predominant component of most connective

tissues within mammalian body. The collagen superfamily comprises 28 known types, type I is the most abundant in the dermis, which constitutes at least 80-90% of the entire collagen [5].

The complex three-dimensional structure of this protein, allows the formation of a scaffold that guarantees a structural and functional support to the dermis and to the overlying epidermis.

Fibroblasts are the main cell population responsible for the production of collagen fibers, which is continuous, and is balanced by the degradation of old collagen by metalloproteinases [6,7]. As we age the collagen fibrils tend to "age" too, due to the formation of covalent bonds (crosslinking) [8] within the structure that destabilize the molecule, making it less elastic and more fragile. This process increases also with the intense exposure to sunlight [9]. The visible results of collagen aging is sagging skin and wrinkles, especially on the face.

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